

Transport by Air - Declaration of Conformity

We hereby confirm that all our power wheelchairs of the Swiss Viva brand, set forth in in the list, below are equipped with maintenance-free and sealed batteries, which are installed in a secure steel box on the wheelchair. These models and batteries meet the requirements set forth in the 2021 – 2022 Edition of the **ICAO Technical Instructions for the Safe Transport of Dangerous Goods by Air (Technical Instructions)** and the 62nd Edition of **the IATA Dangerous Goods Regulations (DGR)**. The batteries can be disconnected from the circuit with a power isolation switch and therefore removing of the batteries is not necessary. Alternatively unplug the bus cable the joystick unit . The power wheelchairs of the Swiss Viva family can all be set into neutral mode by the switch on the top of the drive motors. Details on the physical dimensions and other technical data find at the data table attached.

Battery type: non-spillable, not restricted for transport by air (IATA/ICO provision 67)

Max. Battery measurements:

Length 306 mm

Width 175 mm

Height 190 mm

Weight 23 kg

Liste der Artikel/Liste d'articles/Article List/Lista degli articoli

Produktname Nome prodotto Product name Nome prodotto	Artikelnummer Numéro de référence Article number Codice articolo	Basic UDI-DI	ab Charge/ Serienr.	Klassifizierung Classification Classificazione	Regel règle rule regola
Swiss Viva	004.40.4000	7649997485006	n/a	I	13
Swiss Viva Plus	004.40.4001				
Swiss Viva Junior	004.40.4002				
Swiss Viva Grand	005.40.4000	7649997485013	n/a	I	13
	005.40.4001				
	005.40.4002				
Swiss Viva Grand S	005.40.4010	7649997485020	n/a	I	13
	005.40.4012	7649997485051			
Swiss Viva Grand M	005.40.4005	7649997485037	n/a	I	13
	005.40.4006				
Swiss Viva Junior S	005.40.4011	7649997485044	n/a	I	13

SKS Rehab AG



Daniel Roth, Managing Director

Schwanden, April 15, 2023

HAUPTNIEDERLASSUNG SCHWEIZ

SKS REHAB AG, Im Wyden, CH-8762 Schwanden

sk@sk-rehab.ch T +41 55 647 35 85, F +41 55 647 35 86

www.sk-rehab.ch

SERVICECENTER DEUTSCHLAND

SKS REHAB AG, Hermann-Kolb-Strasse 35b, D-90475 Nürnberg

sk@sk-rehab.com T +49 911 988 12 112, F +49 911 988 12 113

www.sk-rehab.com

Technical Data ¹		Swiss Viva	Swiss Viva junior	Swiss Viva Grand	Swiss Viva Grand S		Swiss Viva Grand M	Swiss Viva Junior S
		Swiss Viva plus			Heckantrieb Rearwheel	Mittelantrieb Midwheel		
Max. Speed	km/h	6/12	6/12	8/12	8/12	8/12	8/12	6/12
Overall length								
incl. legrests	cm	93	90	98	97	110	110	107
w/o legrests, transportation length	cm	66	66	73	73	n.a.	n.a.	93
Width	cm	60	60	63	63	59/65	59	59
Turning radius incl. legrests	cm	79	76	98	98	65	65	65
Seatwidth standard ²	cm	36 - 53	25 - 35	38 - 50	38 - 48	38 - 48	38 - 50	27 - 38
Seatdepth standard ²	cm	35 - 50	35 - 42.5	40 - 54	42.5-58	42.5-58	40 - 54	30 - 40
Seating height								
w/o seat cushion	cm	45	45	43.5	49.5/51	46.5/48	39	47.5
incl. SKS standard seat cushion	cm	48	48	46.5	52.5/54	49.5/51	42	50.5
Tilt in space								
mechanical tunable	°	0 - 17	0 - 17	n.a.	n.a.	n.a.	n.a.	n.a.
powered	°	0 - 27	0 - 27	0 - 40/50	0 - 30	0 - 40/50	0 - 40/50	0 - 30
Backrest recline								
mechanical tunable	°	90 - 102	90 - 102	82 - 135	n.a.	82 - 135	82 - 135	n.a.
powered	°	90 - 150	90 - 150	82 - 153	90 - 150	82 - 153	82 - 153	90 - 133
Seat lift	cm	0 - 30	0 - 30	0 - 30	0 - 30	0 - 30	0 - 30	0 - 30
(Powered) elevating legrests		✓	✓	✓	✓	✓	✓	✓
Tunable legrest carrier angle	°	n.a.	n.a.	0 - 10	n.a.	n.a.	n.a.	n.a.
Front wheel diameter	cm	20/23/26	20/23/26	23/26	23/26	17.5	17.5	17.5
Rear wheel diameter	cm	35	35	35	35	35	35	35
Driving motor	W	2 x 350*	3 x 350*	2 x 350	2 x 350	2 x 400	2 x 400	2 x 400
Batteries	Ah	2 x 12V / 63 Ah (C20)	2 x 12V / 63 Ah (C20)	2 x 12V / 78 Ah (C20)	2 x 12V / 78 Ah (C20)	2 x 12V / 63 Ah (C20)	2 x 12V / 63 Ah (C20)	2 x 12V / 63 Ah (C20)
Loading capacity ²	kg	120	100	180	120	120	120	70
Max. kerb height	cm	5-6	5-6	8	6-8	8	8	8
¹ Subjet to technical changes		*ab 2020	*ab 2020					
² Custom made measures on demand								

HAUPTNIEDERLASSUNG SCHWEIZ

SKS REHAB AG, Im Wyden, CH-8762 Schwanden

sks@sks-rehab.ch T +41 55 647 35 85, F +41 55 647 35 86

www.sks-rehab.ch

SERVICECENTER DEUTSCHLAND

SKS REHAB AG, Hermann-Kolb-Strasse 35b, D-90475 Nürnberg

sks@sks-rehab.com T +49 911 988 12 112, F +49 911 988 12 113

www.sks-rehab.com