

Transport by Air - Declaration of Conformity

We hereby confirm that all our power wheelchairs of the Swiss Viva brand, set forth in in the list, below are equipped with maintenance-free and sealed batteries, which are installed in a secure steel box on the wheelchair. These models and batteries meet the requirements set forth in the 2021 – 2022 Edition of the **ICAO Technical Instructions for the Safe Transport of Dangerous Goods by Air (Technical Instructions)** and the 62nd Edition of the **IATA Dangerous Goods Regulations (DGR)**. The batteries can be disconnected from the circuit with a power isolation switch and therefore removing of the batteries is not necessary. Alternatively unplug the bus cable the joystick unit . The power wheelchairs of the Swiss Viva family can all be set into neutral mode by the switch on the top of the drive motors. Details on the physical dimensions and other technical data find at the data table attached.

Battery type: non-spillable, not restricted for transport by air (IATA/ICO provision 67)

Max. Battery measurements:

Length 306 mm, Width 175 mm, Height 190 mm, Weight 23 kg

Liste der Artikel/Liste d'articles/Article List/Lista degli articoli

Produktname Nome prodotto Product name Nome prodotto	Artikelnummer Numéro de référence Article number Codice articolo	Basic UDI-DI	ab Charge/ Serienr.	Klassifizierung Classification Classificazione	Regel règle rule regola
Swiss Viva Swiss Viva Plus Swiss Viva Junior	004.40.4000 004.40.4001 004.40.4002	7649997485006	n/a	I	13
Swiss Viva Grand	005.40.4000 005.40.4001 005.40.4002	7649997485013	n/a	I	13
Swiss Viva Grand S	005.40.4010 005.40.4012	7649997485020 7649997485051	n/a	I	13
Swiss Viva Grand M	005.40.4005 005.40.4006	7649997485037	n/a	I	13
Swiss Viva Junior S	005.40.4011	7649997485044	n/a	I	13
Swiss Viva Grand F	005.40.4020	7649997485068	n/a	I	13

SKS Rehab AG



Daniel Roth, Managing Director

Schwanden, June 30, 2024

HAUPTNIEDERLASSUNG SCHWEIZ

SKS REHAB AG, Im Wyden, CH-8762 Schwanden
sks@sks-rehab.ch T +41 55 647 35 85, F +41 55 647 35 86
www.sks-rehab.ch

SERVICECENTER DEUTSCHLAND

SKS REHAB AG, Hermann-Kolb-Strasse 35b, D-90475 Nürnberg
sks@sks-rehab.com T +49 911 988 12 112, F +49 911 988 12 113
www.sks-rehab.com